

KEYPOINT® Focus

Compact EMG/NCS/EP System



Dantec Keypoint Focus is the next generation EMG system - delivering the quality and results you expect, at an economical price.

Its new compact design and intuitive, powerful functionality make the Dantec Keypoint Focus EMG system ideal for space conscious professionals in hospitals and private offices.

- NEW** 3 and 6 channel options
- NEW** Increased performance of amplifiers and stimulators
- NEW** Keypoint.NET version 2.10 software now with tremor analysis
- NEW** Microsoft Windows® 7
- NEW** All-in-One EMG system with dedicated controls, integrated stimulators and high quality audio



HANDLE DIFFICULT CASES MORE QUICKLY AND EFFICIENTLY

- Protocol-based testing
- Broad range of powerful analysis tools
- Automatic comparison with normative data
- The only system available with a complete MUP normative database
- Fast, automatic report generation

DESIGN AN EMG SYSTEM THAT'S PERFECT FOR YOU

- Select the configuration that works for you - portable notebook systems to fully featured workstations
- Purchase only the applications you need
- Completely control all measurements, data tables/calculations and displays
- Create customized reports

IT-FRIENDLY FEATURES

- Microsoft SQL database
- Automated study data back-up/archiving functionality
- Complete Active Directory Services
- GDT/HL7 connectivity


The all-new Dantec Keypoint Focus EMG system is based on nearly 50 years of experience in EMG system design. Quality, accuracy, and dependability - what you have come to expect from the Dantec brand.

Versatile EMG/NCS/EP Software

natus
neurology




Multimodality NCS with reference values




Advanced H Reflex



Advanced Decrement test



Multi-MUP with reference values



Peak-triggering SF-EMG



EP with reference values

Please consult www.natus.com for your local sales & service office.

natus

Natus Medical Incorporated
1501 Industrial Road
San Carlos, CA 94070 USA
Tel: +1 650 802 0400
Tel: +1 800 303 0306 (toll-free US/Canada)
Fax: +1 650 802 0401
www.natus.com

natus
europe

Natus Europe GmbH
Bärmannstrasse 38
D-81245 Munich
Germany
Tel: +49 89 83942 0
Fax: +49 89 83942 186
www.natus.com

Manufactured by:
Alpine Biomed ApS
Tonsbakken 16-18
DK-2740 Skovlunde
Denmark
Tel: +45 4457 9000
Fax: +45 4457 9010

Disclosure: See product labeling for specific directions regarding product use, contraindications, warnings, precautions and potential adverse events.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

P/N 007380A
©Natus Medical Incorporated 2011
All Rights Reserved

Keypoint is a registered trademark of Alpine Biomed in the US and other countries. Microsoft and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.